INSTRUCTIONS FOR USE

RECELL GO® Autologous Cell Harvesting Device
AVRL0203, AVRL0204, AVRL0304, AVRL0404, AVRL0504, AVRL0604

The RECELL GO Autologous Cell Harvesting Device (RECELL GO Device) should be used only by licensed healthcare professionals trained in the use of the device.

THIS PAGE INTENTIONALLY LEFT BLANK

A. BACKGROUND	4
A1 DEVICE DESCRIPTION	4
A2 INTENDED PURPOSE	4
A3 INDICATIONS FOR USE	5
A4 INTENDED USERS	5
A5 PATIENT TARGET GROUP	5
A6 CONTRAINDICATIONS	5
A7 WARNINGS	5
A8 PRECAUTIONS	6
A9 CAUTION	6
A10 ADVERSE REACTIONS AND SERIOUS INCIDENTS	7
A11 CLINICAL BENEFITS	7
A12 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE	7
A13 DOSAGE	7
A14 HOW SUPPLIED	8
A15 COMPONENT STERILIZATION AND TESTING	9
A16 DISPOSAL	9
A17 CLEANING	9
B. TREATMENT	10
B1 WOUND BED PREPARATION	10
B2 REQUIREMENTS	10
B3 RECELL GO DEVICE SET-UP	10
B3A REUSABLE PROCESSING DEVICE	11
B3B PREPARE SINGLE USE CARTRIDGE AND COMPONENTS	13
B4 HARVEST SKIN	17
B5 ADD DONOR SKIN TO THE CARTRIDGE	19
B6 PREPARE CELL SUSPENSION	20
B7 DELIVER CELL SUSPENSION	21
B8 CLEANING OF THE PROCESSING DEVICE	23
B9 STORAGE OF THE PROCESSING DEVICE	25
B10 USER INTERFACE	26
C. AFTERCARE	31
C1 SUBSEQUENT DRESSINGS	31
C2 AFTERCARE PRECAUTIONS	31
C3 SCAR MANAGEMENT	32
D. SYSTEM SPECIFICATIONS	32
D1 TECHNICAL SPECIFICATIONS	32
D2 OPERATION AND STORAGE CONDITIONS	32
D3 INTENDED USE ENVIRONMENT	33
D4 ESSENTIAL PERFORMANCE	33
E. ELECTROMAGNETIC COMPATIBILITY	34
F. TROUBLESHOOTING	36
G. APPLICABLE SYMBOLS	37

A. BACKGROUND

A1 DEVICE DESCRIPTION

RECELL GO® Autologous Cell Harvesting Device comprises two separate units. The RECELL GO Processing Device is a durable, powered, reusable device that provides the control for the sterile single-use RECELL GO Cartridge. The sterile, single-use RECELL GO Cartridge is part of the RECELL GO Preparation Kit which consists of components used for preparing the RECELL Spray-On Skin® Cells. The RECELL GO Preparation Kit component boxes include the RECELL Enzyme, Buffer solution, Sterile water, Nozzles, Blunt Fill Needles, 10ml Disposable Syringes and Guarded Disposable Scalpel.

The clinician places split-thickness skin samples into the single-use cartridge along with RECELL Enzyme™ and delivery solutions. The RECELL GO Cartridge locates into the RECELL GO Processing Device which then performs both enzymatic and mechanical disaggregation. Sterile syringes and blunt needles are provided to draw up the RECELL Spray-On Skin Cells from the cartridge once the processing has completed. Sterile nozzles are provided to then apply the RECELL Spray-On Skin Cells to the patient's prepared treatment sites at the point-of-care.

The cell suspension contains a mixed population of cells, including keratinocytes, fibroblasts, and melanocytes, obtained from the disaggregation of the skin sample. Additionally, sub-populations of keratinocytes critical for re-epithelialization have been identified in the cell suspension including basal keratinocytes, suprabasal keratinocytes, and activated keratinocytes. The Enzyme used to process the cells is a biological agent and as such may have slight variations in color and texture.

A2 INTENDED PURPOSE

The RECELL GO Autologous Cell Harvesting Device is intended to be used to disaggregate cells from a patient's split-thickness skin sample and to collect these cells for autologous application to a prepared wound bed.

The RECELL GO Device consists of the following components:

- **RECELL GO Processing Device (RPD):** The RECELL GO Processing Device is a durable, powered, reusable device that is intended to provide the automation for the sterile single-use RECELL GO Cartridge.
- **RECELL GO Processing Cartridge (RPC):** The RECELL GO Cartridge is a sterile, single-use component of the RECELL GO Preparation Kit, which consists of additional components intended to be used for preparing the RECELL Spray-On Skin Cells.
- **RECELL Enzyme:** The Enzyme is intended to provide enzymatic break-down of the extra- cellular matrix of skin sample to prepare cell suspension.
- **Buffer:** The buffer solution is intended to be used as a carrier medium for the cell suspension (temporarily support cell viability).
- **Sterile Water:** The sterile water is intended to reconstitute the Enzyme.
- **Nozzles:** The Spray Nozzles attach to the filled syringe and intended to provide optimal delivery of the cell suspension over the treatment area.

- **Blunt Fill Needle:** The Blunt Fill Needle is intended to draw up fluids from the RECELL Preparation Kit (RPK) component vials and transfer the fluid to the RECELL Processing Cartridge (RPC). The needle is also intended to draw up processed cells from the RPC and transfer the fluid or filtered cell suspension into the delivery syringe.
- **10 ml Disposable Syringe:** The Syringe is intended to transfer the sterile water to reconstitute the Enzyme, transfer the Enzyme and buffer solution to the RPC, and deliver the final cell suspension to the patient.
- **Guarded Disposable Scalpel:** The scalpel is intended to cut skin samples into pieces prior to placing in the RECELL GO Cartridge.

A3 INDICATIONS FOR USE

The cells can be used for autologous application to the prepared wound bed as determined by the physician for a variety of indications such as:

- Burns, or other acute wounds
- Scars
- Chronic wounds or other dermal defects.

A4 INTENDED USERS

Healthcare Professionals trained on the steps involved for safe and effective use of RECELL GO Processing Device. Successful completion of the training program is required prior to allowing the site to use the device without a company representative present.

A5 PATIENT TARGET GROUP

Patients of any age group weighing 3.5 kg and above with need for autologous skin cell suspension application for Burns, or other acute wounds; Scars; Chronic wounds or other dermal defects.

A6 CONTRAINDICATIONS

- RECELL GO is contraindicated for the treatment of wounds clinically diagnosed as infected or with necrotic tissue present in the wound bed.
- RECELL GO is contraindicated for the treatment of patients with a known hypersensitivity to trypsin or compound sodium lactate solution (Hartmann's Solution).
- The skin sample collection procedure specified for use of RECELL GO should not be used with patients having a known hypersensitivity to anesthetics, adrenaline/epinephrine, povidone-iodine, or chlorhexidine solutions.

A7 WARNINGS

- Autologous use only.
- Control infections on wounds prior to application of the cell suspension.
- Excise the necrotic tissues on wound bed prior to application of the cell suspension.
- Wound beds treated with a cytotoxic agent (e.g., silver sulfadiazine) should be rinsed prior to application of cell suspension.
- The RECELL GO Preparation Kit cartridge and components are provided to the healthcare professional sterile and are intended for single use.
- Do not reuse, freeze, or re-sterilize device components.
- Handle using aseptic technique.

- Do not use RECELL GO Preparation Kit or device components if packaging is damaged or there are signs of tampering.
- Choose a skin sample donor site that shows no evidence of surrounding cellulitis or infection.
- For optimum cell viability, the skin sample should be processed immediately after harvesting.
- The Enzyme is derived from animal tissue and, although strict controls have been implemented in the manufacturing process to minimize the risk of pathogen contamination, a small risk of contamination exists and absolute freedom from infectious agents cannot be guaranteed.
- Contaminated materials and waste must be disposed of using appropriate biohazard waste receptacles.
- This device is intended to produce heat as part of its normal operation. Avoid touching the heating plate directly.
- The device is not intended to be calibrated or inspected (as part of preventative maintenance) in the field. Unauthorized device modification may impact the safety of the user.
- The use of RECELL GO Processing Device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Although RECELL GO Processing Device is designed for electromagnetic immunity, use of cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Although RECELL GO Processing Device is designed to be unaffected by typical electrostatic discharge (ESD), very high levels of ESD can result in a temporary suspension of normal operation requiring the operator to recover the operations.

A8 PRECAUTIONS

- RECELL GO is not intended to be used alone (i.e., without meshed autograft) for treatment of acute full-thickness burn wounds or full-thickness skin defects after traumatic avulsion (e.g., degloving) or surgical excision (e.g., necrotizing soft tissue infection) or resection (e.g., skin cancer).
- If the harvested skin sample is processed at ambient temperature above 43°C, it may affect the effectiveness of the enzyme.
- The RECELL GO Processing Device shall be connected to a grounded mains outlet with a supplied hospital-grade power cord.
- Using cleaning agents other than non-woven wipes could damage the equipment or pose an electrical hazard.
- Only use components or accessories provided as part of the RECELL GO Device and follow the manufacturer's instructions for use to ensure optimal device performance.

A9 CAUTION

- The device is intended for use only in a healthcare professional environment complying to electromagnetic limits for medical devices contained in IEC 60601-1-2.
- Follow standard precautions including use of appropriate Personal Protective Equipment (PPE) while cleaning and using the processing device.
- Do not use sterile RECELL GO Preparation Kit components beyond the stated expiration date indicated on the adhesive Lot # and Expiration Date labels on the outer box packaging.
- Do not operate the device while in transit.

A10 ADVERSE REACTIONS AND SERIOUS INCIDENTS

Any adverse reaction or suspected adverse reaction related to RECELL GO should immediately be reported to AVITA Medical, and any serious incident that has occurred in relation to the device should be reported to AVITA Medical and the competent authority of the Member State in which the user and/or patient is established.

Global Customer Service phone number for AVITA Medical: (+1) 833.462.8482.

A11 CLINICAL BENEFITS

RECELL GO is used at a patient's point-of-care to prepare an autologous cell suspension using minimal donor skin. Application of the cellular suspension promotes definitive closure and restoration of natural pigmentation in cutaneous wounds and defects.

A12 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of safety and clinical performance will be available in the European database on medical devices (EUDAMED) public website: <https://ec.europa.eu/tools/eudamed>

A13 DOSAGE

The RECELL GO System consists of the reusable, powered RECELL GO Processing Device and the sterile single-use RECELL GO Preparation Kit consisting of the RECELL GO Cartridge and Disaggregation Head, RECELL Enzyme plus sterile items contained in the component boxes for the RECELL GO Accessories and Procedure Packs.

The contents of each RECELL GO Preparation Kit are sufficient to prepare up to 24 ml of cell suspension which can be used to cover a wound area of up to and including 1,920 cm².

A14 HOW SUPPLIED

The RECELL GO Autologous Cell Harvesting System consists of:


RECELL GO Processing Device (Model #: AVRL0204, AVRL0304, AVRL0404, AVRL0504, AVRL0604)

- 1 x RECELL GO Processing Device
- 1 x Power Cord

RECELL GO Preparation Kit (Model # AVRL0203)

- **RECELL GO Cartridge Tray (UDI-DI #: 00851707004378)**
 - 1 x RECELL GO Cartridge
 - 1 x Disaggregation Head
- **RECELL GO Accessories (UDI-DI #: 00851707004385)**
 - 1 x RECELL Enzyme
 - 1 x 10 ml Buffer
 - 1 x 30 ml Buffer
 - 1 x 10 ml Sterile Water
 - 2 x Nozzle Pouches (4 Spray Nozzles)
- **RECELL GO Procedure Pack A (UDI-DI #: 00851707004392)**
 - 8 x Blunt Fill Needles
 - 4 x 10 ml Syringes
- **RECELL GO Procedure Pack B (UDI-DI #: 00851707004408)**
 - 4 x 10 ml Syringes
 - 1 x Scalpel

RECELL GO System Components

Upon receiving RECELL GO, examine the packaging for external signs of damage. If the external packaging or the packaging for any of the individual components appears damaged, contact your AVITA representative immediately. Do not use any components of the device if the packaging appears damaged. If returning RECELL GO, ensure all original packaging and components are returned with the device.

A15 COMPONENT STERILIZATION AND TESTING

- The RECELL GO Cartridge Tray and Syringes have been sterilized by ethylene oxide.
- The RECELL Enzyme has undergone filtration and terminal sterilization by gamma irradiation.
- The Nozzle Pouches, Blunt Tip Needles, and Scalpel have been sterilized by gamma irradiation.
- The Buffer and Sterile Water have been sterilized using steam.

A16 DISPOSAL

- The RECELL GO Preparation Kit is intended for single use. The RECELL GO Preparation Kit is not reusable and should be discarded after use. Reuse may lead to infection or disease transmission.
- Follow local regulations for proper disposal.
- Contaminated materials and waste must be disposed of using appropriate biohazard receptacles.
- The RECELL GO Processing Device shall be disposed of as per the electrical and electronic equipment disposal process established by the healthcare facility in order to avoid harm to the environment.

A17 CLEANING

It is highly recommended to clean the RECELL GO Processing Device after retrieving from storage and in between uses as part of regular maintenance. Refer to Section B8 for cleaning instructions.

B. TREATMENT

B1 WOUND BED PREPARATION

Clean, vascularized wound bed – To optimize the treatment, the cell suspension should only be applied to a clean, vascularized wound bed with no remaining necrotic tissue. This can be achieved with either dermabrasion using a rotating diamond-head burr, laser ablation, sharp dissection, or other alternative techniques, depending on the nature of the wound.

Infection free – The cell suspension must not be used in the presence of any contamination or infection, as initial re-epithelialization and long-term viability are highly dependent on the absence of infection. Prophylactic antibiotics may be prescribed if the patient is at risk of contamination or infection. Wound swabs for up to-date microbiology are recommended 48 hours prior to the planned surgery.

Pinpoint bleeding – The wound bed should be prepared so that pinpoint bleeding is observed. Accurate debridement to the level of viable tissue is essential; all necrotic tissue must be removed.

B2 REQUIREMENTS

The following materials and instruments will be needed during the procedure:

- Sterile field
- Non-sterile preparation area
- Procedure table
- Personal Protective Equipment
- Skin preparation solution
- Local anesthetic with adrenaline where not contraindicated
- Wound bed preparation tool of choice
- Skin harvesting instrument, e.g., dermatome
- Sterile drape
- Sterile gloves
- Sterile gauze
- Sterile saline
- Sterile blunt forceps
- Sterile ruler, marker
- Appropriate wound dressings - see “Apply Dressings” in Section B7 for details.

B3 RECELL GO DEVICE SET-UP

Perform the following set-up steps in the order shown to avoid setup errors. A procedure guide describing the set-up process is included with the System for reference during a procedure.


Note: Follow software prompts to load RECELL GO Cartridge with the appropriate components.

B3A REUSABLE PROCESSING DEVICE

The RECELL GO Processing Device may be used for (up to) 200 uses with separate single-use RECELL GO Cartridges. The device screen will indicate the number of RECELL GO preparations. Once the count reaches 180 of 200, the RECELL GO Processing Device is within 10% of reaching its limit. It is recommended that you contact your AVITA Medical representative to discuss the appropriate time to replace your RECELL GO Processing Device in order to ensure uninterrupted patient care. It is important to remember that your RECELL GO Processing Device will not prepare any RECELL GO Cartridges beyond the 200-use limit. Once 200 uses have been reached, the RECELL GO Processing Device **MUST** be replaced to maintain effective cell harvesting.

FIRST TIME USE: Remove the RECELL GO Processing Device from the packaging along with the power cable.

Ensure the RECELL GO Processing Device is placed into a clean and flat area and plugged in to a power socket.

CAUTION: The device is intended for use only in a healthcare professional environment complying to electromagnetic limits for medical devices contained in IEC 60601-1-2.

Switch the device ON at the back. Upon a successful test, the screen will now display:

Total RECELL Cartridges Used XX of 200


The system will then ready itself. Please wait and don't interact with the system at this time.

Preparing RECELL Processing Device
Please Wait

The screen will now display:

PLEASE CLEAN

To clean the RECELL GO Processing Device, select the flashing orange CLEAN button.

The screen will now display:

WIPE CLEAN

The shaft will lower and the device is ready for cleaning.

Refer to Section B8 for cleaning instructions.

The screen will display:

Select DONE When Complete

Press the flashing green Done button after the cleaning is completed and the shaft will retract back up.

The screen will now display:


GO PREPARE
RECELL Cartridge With Enzyme, Buffer and Skin For Processing

Move to the next step to start preparation of the sterile RECELL GO Cartridge.


B3B PREPARE SINGLE USE CARTRIDGE AND COMPONENTS

Open the RECELL GO Preparation Kit and place the Procedure Guide and the component boxes containing the RECELL GO Accessories, including the pouched RECELL Enzyme and Procedure Pack items in a non-sterile area. Place the RECELL GO Cartridge Tray in a clean area. The outer shelf box may be discarded.

Remove the lid from the RECELL GO Cartridge Tray and aseptically transfer the RECELL GO Cartridge and Disaggregation Head onto a sterile draped procedure table.


Check the underside of the RECELL GO Cartridge to ensure proper alignment and the pins are centered.


PREPARE ENZYME (WELL A)

- In the non-sterile work area, remove the RECELL Enzyme from the pouch. Remove the cover from the vial marked Enzyme to expose the injection diaphragm. The diaphragm of the Enzyme vial may be wiped with a sterile alcohol wipe and allowed to dry; however, this step is optional.
- Connect a sterile blunt fill needle to a sterile 10 ml syringe and draw up the entire volume of sterile water.
- Inject the entire volume of sterile water into the Enzyme vial. **DO NOT USE** buffer at this stage as this may inhibit the Enzyme action.

- Mix gently until dissolved. Do not shake; use care to avoid foaming. Draw the Enzyme back into the syringe.

- Using aseptic technique, dispense the entire volume of Enzyme into Well A. Discard syringe and needle.

NON-STERILE PREPARATION AREA

STERILE AREA

PREPARE BUFFER (WELL B)

The two (2) buffer vials are to remain in the non-sterile preparation area.

Using aseptic technique introduce the following items into the sterile field.


- 1 x 10 ml Syringe
- 1 x Blunt Fill Needle

- Remove the cover from the 10 ml buffer vial. The diaphragm of the buffer vial may be wiped with a sterile alcohol wipe and allow to dry; however, this step is optional.

- Attach the sterile needle to a clean 10 ml syringe.

- Hold the buffer vial to allow for the volume to be drawn up by person in the sterile field.


HOLD

DRAW UP

- Draw up the entire volume of buffer (approximately 10 ml) from the vial being held by the person in the non-sterile area.

- Dispense the entire volume of buffer into Well B.


NON-STERILE PREPARATION AREA	STERILE AREA
PREPARE BUFFER (WELL C)	
Use the same 10 ml syringe and needle to draw up the relevant volume of buffer from the 30 ml vial (refer to Table 1) and add to Well C of the Cartridge.	

Table 1 - Guidance on Skin Sample and Buffer Volume for Spray-On Skin Cell Preparation

Treatment Area	Skin Sample Size	Starting Volume of Buffer	Approximate Resultant Spray-on Skin Cells Volume
up to 320 cm ²	2 cm x 2 cm (4 cm ²)	4.5 ml	4 ml
up to 480 cm ²	3 cm x 2 cm (6 cm ²)	6.5 ml	6 ml
up to 960 cm ²	2 ea. 3 cm x 2 cm (12 cm ²)	12.5 ml	12 ml
up to 1440 cm ²	3 ea. 3 cm x 2 cm (18 cm ²)	18.5 ml	18 ml
up to 1920 cm ²	4 ea. 3 cm x 2 cm (24 cm ²)	24.5 ml	24 ml

NON-STERILE PREPARATION AREA	STERILE AREA
PREPARE DELIVERY SET ITEMS (WELL D)	
Introduce all of the remaining Accessories and items from Procedure Pack A and B to the sterile field. 6 x 10 ml Syringes 6 x Blunt Fill Needles 4 x Spray Nozzles 1 x Scalpel	

B4 HARVEST SKIN

Skin Sample Type

It is essential that the skin sample harvested is a thin, split-thickness skin sample that leaves pinpoint bleeding at the donor site. The thickness of the skin sample will vary with the body site and patient age and should be in the range of 0.006 – 0.008 in (0.15 – 0.20 mm). The use of a dermatome, or similar device is recommended.


Size of Skin Sample

Choose the appropriate skin sample size for the application. Each square centimeter skin sample can create 1 ml of cell suspension for treatment of an area of up to 80 cm². Each 6 cm² (3 cm x 2 cm) skin sample can yield approximately 6 ml of cell suspension; each RECELL GO Device can process up to four 6 cm² skin samples for a maximum of 24 ml of cell suspension. This can be used to treat an area of approximately up to 1,920 cm².

Table 1 provides guidance for skin sample needed for several treatment area sizes.

Choice of Donor Site

It is essential the donor site is clean, of appropriate depth, and shows no evidence of surrounding inflammation or infection. Choose a donor site of glabrous tissue when creating suspension for glabrous tissue regeneration.

Using the preferred instrument such as a dermatome, take a split-thickness skin sample from the donor site of thickness 0.006 – 0.008 in (or 0.15 – 0.20 mm).

The skin sample may be trimmed from skin harvested for split-thickness skin grafting. Use Table 1 to estimate the skin sample size needed or calculate by taking 1/80 of the total treatment area.

Clean the donor site with antiseptic solution such as povidone-iodine or chlorhexidine. Allow the antiseptic to dry before removing with sterile saline (antiseptic solutions may be cytotoxic and as such, may affect cell viability if left on the skin sample site).

If desired, infiltrate the subcutaneous tissue with a tumescent solution of choice, to provide a firmer surface and anesthesia for taking the skin sample. Ensure that anesthetic is not injected intradermally.


The donor site area may be lubricated (e.g., with sterile mineral oil) to ease travel of the dermatome.

Due to the thick keratin layer found on glabrous skin, it is necessary to take two shaves over the same site in these areas. Discard the first sample and process the second skin sample to create the cell suspension.


Use the scalpel to cut skin samples into pieces no larger than 3 cm x 2 cm prior to placing in the RECELL GO Cartridge. A maximum of 24 cm² can be placed in Well A. Keep the skin samples moist in sterile gauze moistened with sterile saline prior to use.

B5 ADD DONOR SKIN TO THE CARTRIDGE


Place donor skin sample(s) into Well A.


Place donor skin sample(s) separately from each other in Well A. It does not matter which side of the skin is placed up in Well A.

Use **blunt** forceps to spread out the skin along the base of Well A. **Do not use pointed forceps so as to not pierce the filter in Well A.**

Gently place the disaggregation head on top of the skin sample in Well A.

Once the RECELL GO Cartridge is prepared, press the flashing green CONTINUE button.

Next the RECELL GO Processing Device screen will display:

VERIFY
RECELL Cartridge Includes Enzyme, Buffer, Skin and Disaggregation Head

Once you have verified the contents of the RECELL GO Cartridge, press the flashing green CONTINUE button.


B6 PREPARE CELL SUSPENSION

The screen will now display:

INSERT
RECELL Cartridge and Press START

A non-sterile assistant should open the door of the RECELL GO Processing Device.

While maintaining aseptic technique, keep the RECELL GO Cartridge level and with Well A containing the Disaggregation Head at the back of the RECELL GO Cartridge, slide it into the RECELL GO Processing Device until the side tabs on the Cartridge meet the Device.

A non-sterile assistant should close the door and press the flashing green START button to start the skin processing sequence. The screen will display:

RECELL Spray-On Skin Cells Processing
XX Minutes Remaining

The buttons will cease flashing and screen will count down until the RECELL Spray-On Skin Cells are ready for collection.


This will take around 35 minutes. Do not open the door while the RECELL Spray-On Skin Cells are processing.

B7 DELIVER CELL SUSPENSION

Prepare the Dressings

Prior to applying the cell suspension, ensure the dressings are cut and prepared for immediate application. The primary dressing may be fixed using surgical glue, sutures, or staples, or held at the lower aspect of the wound prior to applying the cell suspension to reduce runoff.

Draw up the Prepared Cell Suspension


An audible alert will sound, and the screen will display:

Process Complete
Aseptically Draw Up The RECELL Spray-On Skin Cells From Well D

The RECELL Spray-On Skin Cells are ready to be drawn aseptically from Well D.

Aseptically transfer the remaining blunt fill needles and syringes into the sterile draped procedure field. Attach a blunt fill needle to the 10 ml syringe. Return to the RECELL GO Processing Device noting the device is not sterile. A non-sterile assistant should open the door of the device and move the RECELL GO Cartridge forward on the base of the device so that Well D is accessible using aseptic technique. The RECELL GO Cartridge can also be placed in a clean area next to the RECELL GO Processing Device, if desired.

Using aseptic technique, seat the needle in the dimple at the base of Well D. **Mix the cell suspension by drawing up and returning back to Well D.** Draw up the required volume to apply to the treatment area using the additional syringes and needles in the sterile draped procedure table as required.

Apply Cell Suspension to the Wound Bed

The cell suspension can be sprayed or dripped onto the wound bed, with the technique (i.e., spraying vs. dripping) dependent on the volume of cell suspension to be applied and size of wound bed.

Prior to application, invert the syringe several times to ensure an even suspension.

Spray Application - Application of greater than or equal to 2 ml of cell suspension

The spray application technique should only be used when there is greater than or equal to 2 ml of cell suspension in the syringe.

Remove the needle from the syringe containing the cell suspension. Locate the nozzle in the sterile draped procedure table. Using firm pressure, attach the supplied spray nozzle to the syringe. Be careful not to press the plunger.

Check that the aperture of the attached spray nozzle faces the wound. Hold the spray applicator approximately 10 cm from the most elevated point of the wound and in a position, such that the first drop of suspension falls onto the wound surface. Apply moderate pressure to the plunger of the syringe. Start spraying at the most elevated part of the wound so that any run-off helps to cover the more dependent areas of the wound. A fine mist of cell suspension should be delivered to the wound surface. To cover a larger area, carefully move the spray applicator in one continuous motion from one side of the wound to the other as you spray.

Drip Application - Application of less than 2 ml of cell suspension

The drip application should be used any time that the remaining volume of cell suspension in the syringe is less than 2 ml. Please note that if the remaining volume of cell suspension is less than 2 ml, there is an insufficient amount of cell suspension to cover a wound that is $\geq 160 \text{ cm}^2$.

Do not remove the needle from the syringe containing the cell suspension.

Starting at the most elevated point of the wound, carefully drip the cells onto the wound surface so that any run-off helps to cover the more dependent areas of the wound.

Note: Following application, it is typical to observe run-off of the suspension from the treated wound; results from prospective randomized clinical studies indicate sufficient cellular attachment is obtained post-application of the cellular suspension as epidermal regeneration for definitive closure is achieved.

Apply Dressings

After applying the cell suspension, cover the wound with a non-absorbent, small pore dressing. Always follow the instructions as set by the dressing manufacturer. Dry dressings may be applied moist at the direction of the healthcare professional by lightly soaking the dressing in sterile saline before dressing the wound. The dressing may be fixed to the wound with surgical glue, sutures, or staples, as necessary.

Secondary dressings that are moderately absorbent, minimally adherent, low shear, and readily removable should be placed over the primary dressing followed by absorbent gauze. Use of known cytotoxic medication (for instance, silver sulfadiazine) is contraindicated for areas treated using RECELL GO. Additional absorbent gauze for padding, as well as a crepe or compression bandages, may be used.


B8 CLEANING OF THE PROCESSING DEVICE

Once the procedure is complete, remove the RECELL GO Cartridge from the RPD and dispose of it per the facility protocol. It is highly recommended to clean the equipment after retrieving and using the device from storage and in between treatment use as part of regular maintenance.

After pressing the flashing green CONTINUE button, the screen will display:

PLEASE CLEAN

To clean the RECELL GO Processing Device, select the flashing orange CLEAN button.

The screen will now display:

WIPE CLEAN

The shaft will lower and the device is ready for cleaning.

Wipe the equipment with non-woven towelettes containing surface agents (e.g., disinfecting wipes) that are non-abrasive, non-corrosive and safe for all non-porous surfaces.

Wipe the entire device (interior and exterior) until all gross soil is visually eliminated including, but not limited to, the following areas:

- Inner chamber where the cartridge is inserted
- Shaft which extends from the top of the inner chamber
- LCD screen and buttons
- Door & door hinges

Utilize as many towelettes as required to perform the wipe down and attention shall be paid to areas such as crevices, seams, fluid spillage and stained areas. Ensure all gross soil has been eliminated before proceeding.

Use a dry lint-free cloth to dry the interior and exterior of the device.

The screen will display the following message with flashing green DONE button:

Select DONE When Complete

Press the green Done button after the cleaning is completed and the shaft will retract back up and the screen will return to display status of RECELL GO Cartridges used. The system is now ready to be stored or used with the next patient.

CAUTION: Follow standard precautions including use of appropriate Personal Protective Equipment (PPE) while cleaning and reusing the processing device.

B9 STORAGE OF THE PROCESSING DEVICE






To store the Processing Device:







1. Turn ON/OFF switch to the OFF position.
2. Remove the power cord from the supply mains and detach the power cord from the device.
3. Store the Processing Device within ambient storage conditions.



CAUTION: Do not operate the device while in transit.

B10 USER INTERFACE




Operating Mode Screen Language







Screen Message	Button Function	Purpose	Response
	None	The system is powered on and running a self-test before proceeding.	None required
	CONTINUE button flashes GREEN	Displays total number of RECELL GO Cartridges used to date. Numbers will be highlighted in GREEN if count is less than 180.	Select CONTINUE to proceed to next step.
	CONTINUE button flashes GREEN	Displays total number of RECELL GO Cartridges used to date. Numbers will be highlighted in ORANGE if count is 180 or greater.	Contact your AVITA Medical representative to discuss when to replace your RECELL GO Processing Device. Select CONTINUE to proceed to next step.
	None	The RECELL GO Processing Device is undergoing a homing process to prepare for the next step. Note: This message will appear after pressing "CONTINUE" at Total RECELL Cartridge Used screen or if the device is recovering from an error state.	None required
	CLEAN button flashes ORANGE and SKIP Buton flashes WHITE.	Reminder to CLEAN the RECELL GO Processing Device before proceeding.	Select CLEAN to proceed with cleaning the RECELL GO Processing Device. Select SKIP if cleaning not necessary.






Screen Message	Button Function	Purpose	Response
	None	The RECELL GO Processing Device prepares for you to wipe it clean by lowering the shaft.	Wipe down the RECELL GO Processing Device as instructed in the IFU and Procedure Guide.
	DONE button flashes GREEN	When cleaning is complete, retracts the shaft before proceeding.	Select DONE once you have completed cleaning the RECELL GO Processing Device as instructed.
	CONTINUE button flashes GREEN	Reminder to step away from RECELL GO Processing Device and prepare RECELL GO Cartridge using the RECELL GO Preparation Kit components and skin.	Select CONTINUE once the contents of the RECELL GO Cartridge have been prepared and included.
	CONTINUE button flashes GREEN	Reminder to pause and verify the contents of the RECELL GO Cartridge before placing it in the RECELL GO Processing Device.	Only select CONTINUE once all contents are confirmed to be in the cartridge.
	START button flashes GREEN	Notification that the RECELL GO Processing Device is ready to accept the placement of the RECELL GO Cartridge.	Select START once RECELL GO Cartridge is successfully inserted within the RECELL GO Processing Device chamber.
	None	Notification of the time remaining for Processing of the RECELL Spray-On Skin Cells.	Do not interact with the Device during this process.

Screen Message	Button Function	Purpose	Response
 PROCESS COMPLETE Aseptically draw up the RECELL® Spray-On Skin Cells from Well D. CONTINUE	CONTINUE button flashes GREEN	Notification that the process is complete and the RECELL Spray-On Skin Cells are ready to be aseptically drawn from Well D of the RECELL GO Cartridge.	Follow the instructions in the IFU or Procedure Guide for drawing the suspension from Well D and preparing the syringes for application of the RECELL Spray-On Skin Cells. Select CONTINUE once Cartridge is removed from the Device.
 MAXIMUM REACHED! Contact AVITA Medical Customer Service for a New RECELL® Processing Device	None	Notification that the RECELL GO Processing Device is not able to process any additional Cartridges and must be replaced.	Contact AVITA Medical Customer Service as soon as possible to arrange for a replacement RECELL GO Processing Device.

Error Mode Screen Language

Screen Message	Button Function	Purpose	Response
 ERROR Cartridge NOT Detected. Remove and Reinsert. If Error Persists, Refer to Instructions for Use. ABORT RESUME	ABORT flashes Yellow	Alert that the cartridge is not properly situated or detected by the Device.	Select ABORT if no cartridge or not ready to proceed with Processing.
 ERROR Cartridge NOT Detected. Remove and Reinsert. If Error Persists, Refer to Instructions for Use. ABORT RESUME	RESUME flashes Yellow	Alert that the cartridge is not properly situated or detected by the Device.	Select RESUME if the cartridge is repositioned and detected.
 ERROR Disaggregation Head not found. Remove Cartridge & place Disaggregation Head. Then reinsert Cartridge. ABORT RESUME	ABORT flashes Yellow	Alert that the disaggregation head has not been properly placed in Well A.	Select ABORT if disaggregation head not available.

Screen Message	Button Function	Purpose	Response
	RESUME flashes Yellow	Alert that the disaggregation head has not been properly placed in Well A.	Select RESUME once disaggregation head is placed in Well A.
	ABORT flashes Yellow	Alert when a previously used cartridge is detected. Using a cartridge more than once is not permitted.	Select ABORT if no new cartridge is available.
	RESUME flashes Yellow	Alert when a previously used cartridge is detected. Using a cartridge more than once is not permitted.	Select RESUME if a new cartridge is available and can be used for RECELL Spray-On Skin Cell Processing.
	None	Alert that there is a motor jam or another disruption to the cycle that has halted Processing.	Unplug and turn off the Power to the RECELL GO Processing Device to continue recovery process.
	ABORT flashes Yellow	Alert that the Processing cycle had been interrupted and did not complete.	Select ABORT if Not continuing RECELL GO process.
	RESUME flashes Yellow	Alert that the Processing cycle had been interrupted and did not complete.	Select RESUME to continue RECELL Spray-On Skin Cell Processing.

Screen Message	Button Function	Purpose	Response
	ABORT flashes Yellow	Alert if the cartridge was left in the Device chamber when selecting CLEAN.	Select ABORT if the shaft is jammed on the well housing plate and the cartridge is stuck in place.
	RESUME flashes Yellow	Alert if the cartridge was left in the Device chamber when selecting CLEAN.	Select RESUME after removing the cartridge from the Device chamber.
	Buttons/button are solid Yellow. Button presses are ignored. Continuous audible tone.	Alert that RECELL GO Processing Device is inoperable.	Contact AVITA Technical Support.
	NO button returns system to previous screen.	Confirmation before aborting the process.	Select NO to go back to ABORT/RESUME prompt.
	YES button proceeds with the ABORT request.	Confirmation before aborting the process.	Select YES to ABORT the RECELL GO process.

C. AFTERCARE

The following information, precautions, and notes provide guidelines for care after RECELL GO treatment. Discuss appropriate aftercare with your AVITA representative and provide the patient with aftercare instructions.

C1 SUBSEQUENT DRESSINGS

The outer dressings and compression bandages may need to be changed if exudate levels are high; however, the primary dressing should remain in place for 6-8 days, or as clinically indicated. Protect the primary dressing during secondary dressing changes.

The primary dressing must not be forcibly removed from area(s) to which it is still adhered. Typically, it can be separated (gently peeled back) as new epidermis is formed.

To prevent trauma, any dressing not easily removed should be soaked with an aqueous or oil-based solution prior to removal.

Once the primary dressing has been removed, an appropriate protective dressing should be applied to protect the wound surface.

Do not use dry dressings as protection over blisters or areas of punctate bleeding, as dried exudate could cause newly regenerated epidermis to adhere to the dressing, leading to potential injury upon dressing removal. Instead, use a sterile greasy or paraffin gauze dressing until any blistering or open areas resolve.

Any signs or symptoms of infection or impaired healing at this stage should be recorded and addressed.

C2 AFTERCARE PRECAUTIONS

- Patients should take necessary precautions to prevent the treated area from getting wet while the wound is still open.
- The primary dressing should remain in place as clinically indicated but is typically no longer required after 6-8 days.
- Up to two additional weeks may be needed after initial closure of the treated area for the newly regenerated epidermis to mature and become robust. During this time protective dressings must be worn, particularly on extremities.
- Use of known cytotoxic medication (for instance, silver sulfadiazine) is contraindicated for areas treated using RECELL GO.
- Patients and caregivers should be provided with adequate information and materials for appropriate protection against re-injury during healing and maturation of the treated area.
- Patients should be advised to refrain from strenuous activity.
- Patients should avoid direct sun exposure. A minimum SPF30 and protective clothing should be worn.
- Patients should be counseled about increased risks of skin cancers after thermal burn wounds, and to notify their treating physician of their prior treatment with RECELL GO if they develop skin cancers.

C3 SCAR MANAGEMENT

When the wound has healed, the patient should be advised to continue to protect the area from any surface trauma and to avoid direct sun. Regular use of sunscreen (SPF30) and twice-daily massage with a non-oily skin moisturizer is recommended.

The patient should be advised that the wound area will change over the subsequent weeks and months. The pigmentation and skin texture will continue to mature and improve during this time and the final result may take up to 12 months to be achieved.

Follow-up procedures should follow standard protocols for the specific injury and treatment given.

D. SYSTEM SPECIFICATIONS

This device meets the following standard - IEC 60601-1 edition 3.2 Medical electrical

D1 TECHNICAL SPECIFICATIONS

RECELL GO Processing Device

Mechanical

Dimensions	32.7 cm X 21.5 cm X 34.5 cm (DXWXH)
Weight	7.6 Kg
Case Material	ABS Lustran 348
IP Rating:	IP22

Electrical

Protection	Class I
Rated Power	230 – 264 VA
Supply Frequency	50/60 Hz
Fuse Specification	250 VAC, 2A Slow blow, 100A breaking capacity, 5.2mmX20mm

D2 OPERATION AND STORAGE CONDITIONS

RECELL GO Processing Device

	Operation	Storage
Temperature	15-35°C	Ambient condition
Relative humidity	30 - 60 % (RH)	Ambient condition
Atmospheric pressure	Ambient condition	Ambient condition

Transportation

To move the RECELL GO Processing Device:

1. Turn off the equipment.
2. Disconnect the power cord from the Mains outlet.
3. Lift the equipment from top and bottom to avoid any damage.
4. Use transit cart compliant to healthcare facility regulations for transporting the device.

RECELL GO Preparation Kit

Upon receiving the RECELL GO Preparation Kit, examine the packaging for external signs of damage. If the external packaging or the packaging for the cartridge appears damaged, contact your AVITA representative immediately. Do not use any components of the sterile Cartridge if the packaging appears damaged. If returning RECELL GO, ensure all original packaging and components are returned with the device.

	Storage
Temperature	20-25°C
Relative humidity	Ambient condition
Atmospheric pressure	Ambient condition

The single-use RECELL GO Cartridge and System Components, including the Enzyme, may be stored in a safe and dry environment.

Note: Do not open or use RECELL GO beyond the expiration date listed on the packaging.

D3 INTENDED USE ENVIRONMENT

RECELL GO is intended for use in an operating room of a hospital, surgical center, or surgical suite. However, do not use RECELL GO near active high-frequency surgical equipment, and do not use RECELL GO near RF shielded room of a magnetic resonance imaging equipment where electromagnetic disturbances are high.

D4 ESSENTIAL PERFORMANCE

The degradation of performance or failure of the RECELL GO Processing Device does NOT cause an unacceptable risk to the patient as per the product risk management. The worst-case harm to the patient is limited to reversible skin tissue damage due to the harvesting of skin from a donor site on the patient. Therefore, the device does not support any essential performance functions.

E. ELECTROMAGNETIC COMPATIBILITY

The RECELL GO Processing Device requires special precautions with regard to electromagnetic compatibility (EMC) and should be used in accordance with the information provided in this manual.

The device complies with the requirements of the EMC international standard IEC 60601-1-2 when used in accordance with the electromagnetic environment specified in the following tables.

Guidance and manufacturer's declaration – electromagnetic emissions

Emission Test	Compliance	Electromagnetic environment – guidance
Radiofrequency (RF) emissions CISPR 11	Group 1	The RECELL GO Processing Device does not generate any RF energy for their internal functioning. Therefore, the RF emissions produced by the device shall be very low.
RF emissions CISPR 11	Class A	The RECELL GO Processing Device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	The RECELL GO Processing Device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	N/A

Guidance and manufacturer's declaration – electromagnetic Immunity

Immunity Test Standard	Compliance Level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electromagnetic compatibility (EMC) IEC 61000-4-3	3 V/m	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device. Otherwise, degradation of the performance of this equipment could result.
Power Frequency (50/60 Hz) Magnetic Field, IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity Test Standard	Compliance Level	Electromagnetic environment – guidance
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Electrical fast transient/burst IEC 61000-4-4	± 1 kV for input/output lines ± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, shorts, interruptions, and voltage variations on power supply input lines IEC 61000-4-11	<5% UT [>95% dip in UT] for 0.5 cycle (<5% UT [>95% dip in UT] for 0.5 cycle) 40% UT [60% dip in UT] for 5 cycles (40% UT [60% dip in UT] for 5 cycles) 70% UT [30% dip in UT] for 25 cycles (70% UT [30% dip in UT] for 25 cycles) <5% UT [>95% dip in UT] for 5 s (<5% UT [>95% dip in UT] for 5 s)	Mains power quality should be that of a typical commercial or hospital environment. If the user requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	RF shielded environment including filtering of all cable passing through the shielding, with a minimum shielding effectiveness and filter attenuation of 20 dB.

WARNINGS:

Use of RECELL GO Processing Device adjacent to or stacked with other equipment.

Use of RECELL GO Processing Device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Use of accessories, transducers, or cables not specified.

Although RECELL GO Processing Device is designed for electromagnetic immunity, use of cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Electrostatic Discharge Warning

Although RECELL GO Processing Device is designed to be unaffected by typical electrostatic discharge (ESD), very high levels of ESD can result in a temporary suspension of normal operation requiring the operator to recover the operations.

F. TROUBLESHOOTING

Cartridge not detected

In an event of the RECELL GO Cartridge not making electrical connection to the RECELL GO Processing Device after multiple attempts the following error message shows up on the LCD screen – ‘Cartridge Not Detected. Remove and reinsert. If Error Persists, Refer to Instructions for Use’. In this scenario, clean the electrical pads with non-woven towelettes before the next attempt. If the issue persists, use a new RECELL Go Preparation Kit, and set up the RECELL GO cartridge according to the instructions in section B3B. Using aseptic technique, retrieve the skin sample from Well A of the existing RECELL GO cartridge and transfer to the replacement cartridge.

Enzyme powder does not dissolve completely

Make sure that the Enzyme is mixed well with the sterile water by gently inverting the vial several times. Often a small amount of particulate matter remains undissolved in the reconstituted solution. This does not reduce the activity of the Enzyme.

Do not use buffer to dissolve the Enzyme as it may interfere with the Enzyme action.

Skin sample is too large, or too thin

Take particular care when harvesting the skin sample. It should be a thin (0.006 – 0.008 in, 0.15 – 0.20 mm) split-thickness graft with just a very thin section of dermis. The skin sample of the appropriate thickness will ensure successful disaggregation of cells. The maximum size of each skin sample recommended for use with the RECELL GO Device is 3 cm by 2 cm.

If the skin sample is too large (greater than the maximum recommended), cut it into a smaller size and discard the excess.

If the skin sample is too thin, you should take another skin sample from a DIFFERENT donor site.

Buffer added to Enzyme vial

If buffer, instead of sterile water, is mistakenly added to the Enzyme vial, the Enzyme activity may be inhibited. If buffer is mixed with the Enzyme powder, the Enzyme should be discarded and a new RECELL GO Preparation Kit used.

Nozzle blocked

If the cell suspension is not easily sprayed, the cell suspension may be dripped onto the wound bed. If the cell suspension does not come out at all, the nozzle attached to the syringe may be blocked by unfiltered particles. Filter the suspension and place in a new 10 ml syringe prior to attaching a new spray nozzle.

Insufficient treatment area coverage

If cell suspension is lost in the application process and sufficient coverage of the treatment area was not achieved, take another skin sample, and repeat the process with a new RECELL GO Preparation Kit to create additional cell suspension and complete the treatment.

G. APPLICABLE SYMBOLS

The symbols below and harmonized symbols may be found on the product or product label.

Symbol	Description	Symbol	Description
	Consult instructions for use or consult electronic instructions for use		Refer to instruction manual/ booklet
	Manufacturer		Do not use if package is damaged and consult the instructions for use
	Date of manufacture		Use-by Date (i.e., Expiration Date)
	Do not re-sterilize		Do not re-use
	Sterile		Sterilized using irradiation
	Sterilized using ethylene oxide		Sterilized using steam or dry heat
	Single sterile barrier system		Single sterile barrier system with protective packaging outside
	Double sterile barrier system		Contains biological material of animal origin
	Temperature limit (Operating/ Storage/ Transport)		Humidity limitation (Operating/ Transport)
	Serial Number		Unique device identifier
	Batch Code/ Lot Number		Catalogue number
	Medical device		Importer
	Authorized representative in the European Union		Authorized representative in Switzerland
	European Mark of Conformity		Prescription only
	Evaluated and classified by Eurofins MET Laboratories, a Nationally Recognized Testing Laboratory (NRTL), for compliance with specific U.S. safety standards		Waste Electrical and Electronic Equipment (WEEE)
IP22	Protected against solid foreign objects and vertically falling water drops (at 15°h)		

For further information regarding RECELL GO Autologous Cell Harvesting Device, contact your local sales representative.


AVITA Medical Americas, LLC (Legal Manufacturer)
28159 Avenue Stanford, Suite 220
Valencia, CA 91355-2271, USA
Tel: +1 833 462 8482
Fax: +1 661 367 9180
Email: customerservice@avitamedical.com

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
Tel: +31.70.345.8570
Fax: +31.70.346.7299

MedEnvoy Global B V
Prinses Margrietplantsoen 33 - Suite 123
2595 AM The Hague
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Australian Sponsor
Emergo Australia
Level 20 Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

©2025 AVITA Medical.

AVITA Medical®, the AVITA Medical logo, RECELL®, RECELL Enzyme™, RECELL GO® and Spray-On Skin® Cells are trademarks of AVITA Medical. All other trademarks are the properties of their respective owners.

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AVITA Medical Americas, LLC (produttore legale)